


Millipore®

Filtration, Separation
& Preparation

MERCK

Just SNAP and go!

SNAP i.d.® 2.0 system for
Western blotting and IHC



The life science business of Merck
operates as MilliporeSigma in the
U.S. and Canada.



Multiple slides, multiple blots, multiple conditions.

There's so much room for experimental variability in traditional immunodetection workflows. For your peace of mind – and ours – we designed the SNAP i.d.[®] 2.0 system to streamline your Western blot and immunohistochemistry experiments.

The concept is simple: a vacuum-driven flow of blocking, antibody and washing solutions reduces slide and membrane handling. That means a lot less shaking, dipping, pouring and waiting.

And now you can process multiple blots and slides in parallel, so it's easy to apply consistent conditions across experiments.

SNAP i.d.[®] 2.0 Protein Detection System for Western Blotting

Unlike conventional Western blotting, where diffusion is the primary means of reagent transport, the **SNAP i.d.[®] 2.0 system** applies a vacuum to actively drive reagents through the membrane. This advanced technology promotes antigen binding and thorough washing, enabling you to better optimize your Western blotting conditions.

Key Benefits

- Faster results
- Faster testing of different antibodies
- Higher throughput of Western blots each day




How does the SNAP i.d.[®] 2.0 system work?

The vacuum-driven SNAP i.d.[®] 2.0 system fully exploits three-dimensional reagent distribution and decreases the immunodetection time from hours to minutes using the following mechanisms:


1. The system increases local antibody concentrations at binding sites by using vacuum filtration, driving the antibody-antigen binding reaction forward and shortening incubation times.
2. Vacuum pulls any residual, unbound antibody out of the membrane, lowering background signal.

Advantages of the SNAP i.d.[®] 2.0 system's vacuum transport feature


- Draws reagents through blotting membrane
- Minimizes over-blocking
- Thoroughly flushes membranes instead of just rinsing
- Reduces incubation times



SNAP i.d.[®] 2.0 system in the Western blotting workflow



Faster Blots, Better Signals.
Comparison of the traditional Western blotting protocol relative to SNAP i.d.[®] 2.0 system's 30-minute protocol.



Blot Holder Specifications

	Midi	Mini	Each Half of MultiBlot
Dimensions (cm)	17.8 x 10.2	12.7 x 9.1	11.4 x 6.4
Max Blot Size (cm)	8.5 x 13.5	7.5 x 8.4	4.5 x 8.5

Equivalent performance, a fraction of the time

You frequently need to optimize antibody conditions using multiple antibody dilutions, across multiple sample types and matrices. With the SNAP i.d.[®] 2.0 system, you can

turn around your optimization experiments in a fraction of the time it takes for traditional Western blot detection.

Western blot optimization of anti-Tau-1 antibody in Alzheimer's disease brain and healthy brain samples

Human brain samples from an Alzheimer's disease patient and from a healthy donor were lysed in CytoBuster™ Protein Extraction Reagent (Cat. No. 71009). Samples were serially diluted and separated by SDS-PAGE. Gels were transferred to Immobilon®-P membrane. Blots were processed in the SNAP i.d.[®] 2.0 system using MultiBlot, Mini and Midi frames with their corresponding blot holders. A control blot was processed by standard immunodetection. All blots were blocked with 0.5% NFDM and probed with primary anti-Tau1 (Cat. No. MAB3420) and secondary HRP-conjugated goat anti-mouse (Cat. No. AP124P) using the conditions indicated above. Blots were incubated with Luminata™ Forte HRP substrate and exposed to X-ray film for 15 minutes.

	SNAP i.d. [®] 2.0 MultiBlot		SNAP i.d. [®] 2.0 Mini Blot		SNAP i.d. [®] 2.0 Midi Blot		Standard Western blot detection	
MW	Alzheimer's brain	Healthy brain	Alzheimer's brain	Healthy brain	Alzheimer's brain	Healthy brain	Alzheimer's brain	Healthy brain
188	1 2 3 4 5	1 2 3 4 5	1 2 3 4 5	1 2 3 4 5	1 2 3 4 5 6 7 8 9	1 2 3 4 5 6 7 8 9	1 2 3 4 5	1 2 3 4 5
98								
62								
49								
38								
28								
14								
6								
3								
Blocking	0.5% NFDM for 2 sec		0.5% NFDM for 20 sec		0.5% NFDM for 20 sec		0.5% NFDM for 1 hr	
Primary Antibody	Anti-Tau1 1:1,000 for 10 min		Anti-Tau1 1:1,000 for 10 min		Anti-Tau1 1:1,000 for 10 min		Anti-Tau1 1:5,000 for 1 hr	
Secondary Antibody	Goat anti-Mouse 1:10,000 for 10 min		Goat anti-Mouse 1:10,000 for 10 min		Goat anti-Mouse 1:10,000 for 10 min		Goat anti-Mouse 1:50,000 for 1 hr	
Total Time	< 30 min		< 30 min		< 30 min		3 hr 30 min	

MultiBlot, Mini Blot and Standard Western blot detection

Lane	Concentration (µg)
1	20
2	10
3	5
4	2.5
5	1.25

Midi Blot

Lane	Concentration (µg)
1	20
2	10
3	5
4	2.5
5	1.25
6	0.63
7	0.31, 0.16 and 0.08
8	0.16
9	0.08



Ordering Information

Product Description	Qty/Pk	Cat. No.
Base System		
The SNAP® i.d. 2.0 systems contain everything you need to get started, including the detection base, 2 blot holding frames, 2 blot holders, 2 antibody collection trays, a blot roller and rolling pad, 2 wetting trays, vacuum tubing and a Quick Start User Guide.		
SNAP® i.d. 2.0 System - Mini (7.5 x 8.4 cm)	1	SNAP2MINI
SNAP® i.d. 2.0 System - Midi (8.5 x 13.5 cm)	1	SNAP2MIDI
SNAP® i.d. 2.0 System - MultiBlot (4.5 x 8.4 cm)	1	SNAP2MB3
SNAP® i.d. 2.0 System - Mini and Midi (7.5 cm and 8.5 x 13.5 cm)	1	SNAP2MM
SNAP® i.d. 2.0 System - Mini and MultiBlot	1	SNAP2MB1
SNAP® i.d. 2.0 System - Midi and MultiBlot	1	SNAP2MB2
Components for Western Blotting Procedures		
SNAP i.d.® 2.0 MultiBlot Holding Frame	1	SNAP2FRMB01
SNAP i.d.® 2.0 Mini Blot Holding Frame (single pack)	1	SNAP2FRMN01
SNAP i.d.® 2.0 Mini Blot Holding Frames (double pack)	1	SNAP2FRMN02
SNAP i.d.® 2.0 Midi Blot Holding Frame (single pack)	1	SNAP2FRMD01
SNAP i.d.® 2.0 Midi Blot Holding Frames (double pack)	1	SNAP2FRMD02
SNAP i.d.® 2.0 MultiBlot Holders (includes 2 well blanks)	50	SNAP2BHMB050
SNAP i.d.® 2.0 Mini Blot Holders	100	SNAP2BHMN0100
SNAP i.d.® 2.0 Midi Blot Holders	100	SNAP2BHMD0100
SNAP i.d.® 2.0 Antibody Collection Tray	20	SNAPABTR
SNAP i.d.® Blot Roller	1	SNAP2RL
Blotting Membranes		
Immobilon®-P PVDF, 0.45 µm, 26.5 x 3.75 cm roll	1	IPVH00010
Immobilon®-P PVDF, 0.45 µm, 7 x 8.4 cm sheet	50	IPVH07850
Immobilon®-P PVDF, 0.45 µm, 8.5 x 13.5 cm sheet	10	IPVH08130
Immobilon®-FL PVDF, 0.45 µm, 26.5 x 3.75 cm roll	1	IPFL00010
Immobilon®-FL PVDF, 0.45 µm, 7 x 8.4 cm sheet	10	IPFL07810
Immobilon®-PSQ PVDF, 0.2 µm, 7 x 8.4 cm sheet	50	ISEQ07850
Immobilon®-P PVDF Sandwich, 0.45 µm, 7 x 8.4 cm	20	IPSN07852
Immobilon®-P PVDF Sandwich, 0.45 µm, 8.5 x 13.5 cm	20	IPSN08132
Reagents for Western Blotting		
Luminata™ Forte Western HRP Substrate	500 mL	WBLUF0500
Luminata™ Crescendo Western HRP Substrate	500 mL	WBLUR0500
Luminata™ Classico Western HRP Substrate	500 mL	WBLUC0500
Immobilon® Western HRP Substrate	500 mL	WBKLS0500
Calbiochem® SignalBoost™ Immunoreaction Enhancer Kit	1 kit	407207
Re-Blot Plus Strong Antibody Stripping Solution, 10X	50 mL	2504
bløk®-CH Buffer	500 mL	WBAVDCH01
bløk®-FL Buffer	500 mL	WBAVDFL01
bløk®-PO Buffer	500 mL	WBAVDP001

SNAP i.d.® 2.0 IHC System


SNAP i.d.® 2.0 Protein Detection System for Immunohistochemistry (IHC) introduces a new capability to the innovative, vacuum-driven SNAP i.d.® 2.0 system. The IHC frame and slide holders allow you to block, probe and stain up to 12 tissue slides per frame. Reduced handling time and multiple-slide processing make this system ideal for antibody and protocol optimization.



Key Benefits

- Eliminates the need for pap pens
- Antibodies can be collected and reused
- Slide handling time is significantly decreased
- Less time spent on wash steps
- Parallel processing of multiple slides

Immunohistochemistry workflow includes blocking, antibody incubations, labeling and wash steps, all of which can be streamlined using the SNAP i.d.® 2.0 Protein Detection System for IHC.



Key Features

- Flexibility of multiple slide configurations enables the processing of 1 to 24 slides at a time
- Compatible with standard IHC slides and protocols
- Compatible with diverse tissue preparations, including formalin-fixed or fresh frozen samples
- Test tracker feature on frame cover helps keep track of IHC steps
- Intuitive format
 - Incorporates blocking, washing and antibody incubation and labeling steps
 - Systematizes handling multiple slides without the cost of automation


How does the SNAP i.d.® 2.0 Protein Detection System for IHC work?

With two individually controlled sides, the system base allows for independent, vacuum-driven processing of either one or two IHC frames. Each of the IHC frames can process between 1 to 12 glass slides through independent vacuum ports.

Each slide holder has an injection/recovery port that enables the manual addition, as well as the removal and recovery, of small volumes of antibodies or reagents; reagents can also be flushed using the vacuum feature if conservation is not a priority.


Comparable performance to traditional methods, even in archival tissue

The SNAP i.d.[®] 2.0 IHC system produces comparable staining to traditional protocols, even in archival tissue. In the first example below, the system was used to detect Aquaporin 1 in archival human kidney tissue. Note the characteristic differential staining of kidney proximal tubule epithelium and renal corpuscle. The second example shows NeuN signal in human brain, localized as expected to neuronal nuclei. Despite processing 12 slides in parallel and shortening the handling time and protocol, the staining is robust and consistent, with no blotchy artifacts that sometimes plague autostainers. There is no apparent tissue degradation, as compared with traditional protocols. Classic histological stains, such as hematoxylin, can be applied using the same system.




Detection of Aquaporin 1 in human kidney tissue (formalin-fixed and paraffin-embedded [FFPE]): SNAP i.d.[®] 2.0 IHC System (sections 1–13) vs. Standard IHC protocol (section M1). Fifteen tissue sections (5 μ m) were assembled on FisherBiotech[®] ProbeOn Plus[™] slides. Slides were rehydrated and antigen retrieved (heat-induced epitope retrieval, HIER) using Reveal Decloaker (Biocare Medical, LLC) in a pressure cooker for 15 minutes at 110°C. Thirteen slides were then processed using the SNAP i.d.[®] IHC system, and one was processed using the manual protocol. Blocking was performed by incubating 10 min with Punisher[™] reagent (Biocare). After three washes with TBST, slides were incubated 2 h with Anti-Aquaporin 1 (Cat. No. AB2219, 1:2,000). After three more TBST washes, slides were incubated 10 min with Anti-Rabbit Secondary Antibody (Biocare). Signal was detected using a HRP-DAB detection kit (Biocare). Slides were washed three times with TBST and were counterstained with hematoxylin for 1 min. After three final washes, slides were dehydrated with four 5-minute changes of 100% ethanol, cleared with three changes of xylene and coverslipped with Ecomount medium (Biocare).

SNAP i.d.[®] IHC System



Manual IHC



Detection of NeuN in human cerebral cortex (FFPE): SNAP i.d.[®] 2.0 IHC System (left) vs. standard manual IHC protocol (right). Anti-NeuN (Cat. No. [MAB377](#), 1:2,000) was used to stain sections of human cerebral cortex using the protocols described in the previous figure.

Ordering Information

SNAP i.d.® 2.0 system

Base System

The SNAP® i.d. 2.0 systems for IHC contain everything you need to get started, including the detection base, IHC frame and incubation cover, slide holders, an assembly fixture, vacuum tubing and a Quick Start User Guide.

Product Description	Cat. No.
SNAP i.d.® 2.0 Protein Detection System – Single IHC	SNAP2IHC
SNAP i.d.® 2.0 Protein Detection System – Double IHC	SNAP2IHC2

SNAP i.d.® 2.0 consumables

Product Description	Qty	Cat. No.
SNAP i.d.® 2.0 IHC Frame	1 EA	SNAP2FRIHC
SNAP i.d.® 2.0 IHC Slide Holders	24/pk	SNAP2SH


Reli(Ab)le™ Antibodies for IHC

Reliable results depend on reliable reagents. When you work with reliable antibodies, you can work efficiently and effectively. We guarantee performance because we manufacture performance. We are committed to manufacturing antibodies that perform as specified and as anticipated. Our quality starts at inception and carries through manufacturing, production and distribution, into the lab and onto the bench. From start to finish, we build a platform for reliable performance that provides you with the ultimate confidence in your findings.

Popular IHC Antibodies

Product Description	Qty	Cat. No.
Anti-Actin Antibody, clone C4	100 µL	MAB1501
Anti-APP A4 Antibody, a.a. 66-81 of APP {NT}, clone 22C11	50 µg	MAB348
Anti-Choline Acetyltransferase Antibody	500 µL	AB144P
Anti-GAD67 Antibody, clone 1G10.2	100 µg	MAB5406
Anti-Microtubule-Associated Protein 2 (MAP2) Antibody	100 µL	AB5622
Anti-NeuN Antibody, clone A60	500 µg	MAB377
Anti-NG2 Chondroitin Sulfate Proteoglycan Antibody	100 µg	AB5320
Anti-Olig-2 Antibody	100 µL	AB9610
Anti-Sox9 Antibody	100 µg	AB5535
Anti-Tyrosine Hydroxylase Antibody	100 µL	AB152
Anti-E2A antibody produced in rabbit		SAB4502929
Anti-TWIST1		SAB2108515
ANTI-OCT3/4 (C279) antibody produced in rabbit		SAB1306212
Anti-β-Catenin antibody produced in rabbit		C2206
Anti-SNAI1		SAB2108482
Anti-TGF β1, C-Terminal antibody produced in rabbit		SAB4502954
Anti-SMAD3 antibody produced in rabbit		SAB2702052
Monoclonal Anti-EGF antibody produced in mouse		SAB5300488
Anti-BMP4 antibody produced in rabbit		SAB4200556
Anti-Hepatocyte Growth Factor antibody produced in goat		H7157

Order now and put the most reliable antibodies to work for you.



Copyright © 2017 Merck KGaA. All Rights Reserved. Merck, Milli-Q, Millipore and the Vibrant M are trademarks of Merck. Sigma-Aldrich, Supelco, SAFC and BioReliance are trademarks of Sigma-Aldrich Co. LLC. or its affiliates. SNAP i.d., Immobilon, CytoBuster, Luminata, Calbiochem, Reli(Ab)le and blök are trademarks of Merck. All other trademarks are the property of their respective owners.
merckgroup.com/life-science

PB6745ENEU_DP Ver. 2.0
2018 - 18522
12/2018

© 2018 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

 **fisher scientific**
part of Thermo Fisher Scientific